Updated *eDOCSNL NL CPP Sidebar



User Guide

Based on guidance from the eDOCSNL Program's Clinical Advisory Committee, an updated version of the provincial NL CPP Sidebar, called **"*eDOCSNL NL CPP Sidebar**" is now available in all EMRs:

The improved value of the updated Sidebar is in:

- the addition of new widgets (boxes), such as "INRs (Most Recent 5)".
- the added details that will now display within widgets, such as Goal details.
- the reordering of the widgets; and
- the ability to individually manage selected population-based primary care Clinical Decision Support (CDS) screening triggers/manoeuvres through the use of applying and cancelling Goals. (Please refer to the User Guide entitled *"eDOCSNL Using Cancelled Goals to Disable CDSs for Specific Patients"* for further details.)

For assistance, please reach out to the eDOCSNL Program at info@eodcsnl.ca

		*eDOCSNL NL CP	P Sideb Prin	t Select S	idebar View	as Timeline
Profile - Curre	nt and Pas	t Only				Ξ
Medical Curren Surgical Hx Past	t Back Pain Lobotomy -					_
Meds						Ξ
Active AMOXICILLIN 500 r ferrous fumarate 30	ng CAPSULE 0 mg (100 mg	iron) TABLET	01May2 01May2	3-04May 3-30Jul2	23 Oral 3 Oral	
Goals						Ξ
Active Goals (6)	Cancelled G	oals				
Active Goals Goal Name Influenza	Target Last	Value	Last Date 01Jan2023	Next Di 01Jan2	ue 024	
Influenza (Flu)			01Jan2023	01Jan2	024 🔒	
Influenza Vaccine			01Jan2023	01Jan2	024 🔒	
FECAL IMMUNO			01Jan2023	01Jan2	025 🔒	
Mammogram			01Jan2023	01Jan2	025 🔒	
PAP SMEAR			01Jan2023	01Jan2	026 🔒	
Clinical Metric	s					=
CHADS2 Score		6			2023-05	-01
Framingham CVD Score	10-year Risk	7.3			2023-05	-01
Clinical Frailty Sco	re	5. Living V Frailty	Vith Mild		2023-05	-01
PHQ-9 Question 1	Score	9			2023-05	-01
GAD-7 Question 1		7		0/	2023-05	-01
SpO2 Height		90		70 CM	2023-05	-01
BP systolic		110		mm Ha	2023-05	-01
Weight		68.0		kg	2023-05	-01
Glucose (meter)		5.6			2023-05	-01
Pulse		85		bpm	2023-05	-01
BP diastolic		85		Hg	2023-05	-01
vvaist circumferend BMI	e	28.3		cm kg/m^2	2023-05	-01
HbA1C		7.0		Agrif ¹ 2	2021-04 Final	-27
HbA1c/TOTAL HEI	MOGLOBIN	80		%	2017-07 Final	-17



TO FIND AND FAVOURITE THE UPDATED SIDEBAR

To access and Favourite the updated *eDOCSNL NL CPP Sidebar, hover and right-click on any existing Sidebar icon from within a patient chart, hover on "Open Sidebar" and select (left-click) on "Open *eDOCSNL NL CPP Sidebar":



Once selected, the *eDOCSNL NL CPP Sidebar can be Favourited by left clicking on the gray heart:



The *eDOCSNL NL CPP Sidebar can also be set as a User Default Preference in their Profile:

Inchitest.med-access.net/ui/stat	ic/b3fce1e66006e378d5384c5a29d3a1fe.html#/admin/user/detail User Settings: bblake General Groups Passwords Macros Signing Options Preferences Filters Preferred Mgmt History -	Templates Workflows Features Permissions • Contra Preferences
My Profile USER ADMINISTRATION User List User Groups Resources Roles	User Interface Preferences Default Tab: Demographics Chart Summary Full Chart Summary Template Sidebar: ************************************	Scheduling Preferences Schedulable Appt Interval: 10 v mins Default Appt Type

Or, if able to access Template Admin, can be set as a Favourite for a User, specific Group or all Users in the EMR:

Template M	anagement														List New Export Hel
Demog	Visits	Tasks	Meds	Profile	Labs	Invest	Consults	Imm	Goals	Appt	CDS	Filter	Workflow	Dashboards	Obs
Dashboards .															New
Name	Type Sidebar	-												User	Group Site/EMR
Name *eDOCSNL NL C	CPP Sidebar				Des Upd	cription ated sidebar from M	NL CPP (2023)			Type Sidebar	User		B,	2	

TIPS:

• The width of a Sidebar can be changed by the user by dragging the divider line:



• Toggle the Sidebar on and off within Tasks by clicking on the Sidebar icon:



• Make a copy of a Sidebar by hovering, right-clicking and selecting "Copy Sidebar". An icon will be added and the copy will be renamed with the user's name added. The copied Sidebar can be Favourited, renamed and edited as required:





DETAILS OF WIDGETS

Profile – Current and Past Only

Deleted Profile items will no longer be displayed in the Sidebar widget, but can always be accessed via the Profile Tab. You will notice that the Current profile items have been separated from the Past on separate tabs in the widget for ease of view and to prevent the widget from getting too cluttered.



Meds

Displays "Active and Recently Active" (only) and has ROUTE added to the display:

Meds				Ξ
Active				
AMOXICILLIN 500 mg CAPSULE	1 Caps	12Apr23-15Apr23	Oral	

Goals

Multiple new features added:

• Active Goals will display the Target Value, Last Value (if numeric), the Last Date of the result, Next Due date based on interval set (or Unknown or Overdue message) and hovering on the Note icon will display any Reason added to that Goal:

Goals					
Active Goals (5)	Cancelled Goals (2)				
Active Goals Goal Name HbA1c/TOTAL HE	Target Last Value < 7 80	Last Date 17Jul2017	Next Due Overdue		
Influenza (Flu)		01Jan2023	01Jan2024	æ	
FECAL IMMUNOAS	S	01Jan2023	01Jan2025		
Mammogram		01Jan2023	01J	ancer	Screening
PAP SMEAR		01Jan2023	01Janzozo	Ð	Ŀ

• Cancelled Goals will display the Last Value (if numeric), the Last Date of the result and the Cancel Reason:

Goals				Ξ
Active Goals (5)	Cancelled Goals (2)			
Cancelled Goal	S			
Goal Name	Last Value	Last Date	Cancel Reason	
Influenza (Flu)		01Jan2023	No longer indicated	
Mammogram		01Jan2023	Patient declined	

Note: Providers now have the ability to individually manage selected population-based primary care Clinical Decision Support (CDS) screening triggers/manoeuvres through the use of applying and cancelling Goals. (Please refer to the User Guide entitled *"eDOCSNL Using Cancelled Goals to Disable CDSs for Specific Patients"* for further details).

Clinical Metrics

This widget has been renamed (previously Patient Vitals) and now displays the most recent value ("last") and date of 16 key clinical metrics, if they are available in the patient's chart. They are sorted in the order of most recent to oldest. (Note that due to a known technical issue within the software, that the BP systolic and BP diastolic readings are still displayed separately):

Clinical Metrics			Ξ
CHADS2 Score	6		2023-05-01
Framingham CVD 10-year Risk Score	7.3		2023-05-01
Clinical Frailty Score	5. Living With Mild Frailty		2023-05-01
PHQ-9 Question 1 Score	9		2023-05-01
GAD-7 Question 1	7		2023-05-01
Waist circumference	100	cm	2023-05-01
BMI	28.3	kg/m^2	2023-05-01
SpO2	98	%	2023-05-01
Height	154.9	cm	2023-05-01
BP systolic	110	mm Hg	2023-05-01
Glucose (meter)	5.6		2023-05-01
Weight	68.0	kg	2023-05-01
Pulse	85	bpm	2023-05-01
BP diastolic	85	mm Hg	2023-05-01
HbA1C	7.0		2021-04-27 Final
HbA1c/TOTAL HEMOGLOBIN	80	%	2017-07-17 Final

Allergies

Includes more detail on the Allergen, Reaction and Severity:

Allergies
DA Sulfa (Sulfonamide Antibiotics):sulfamethoxazole/trimethop Anaphylaxis Severe
No review: No non-drug allergies.
No review: No drug intolerances.
No review: No non-drug intolerances.

INR Results (Most Recent 5)

A new widget which will display the more recent 5 INR results (if applicable to the patient), including those resulted from Point of Care testing equipment:

INR Results (Most Recent 5)							
Name INR	Value	Units	Date 2023-04-04 Final	Flag			
INR	1.10		2023-03-09 Final				
INR;POINT OF CARE	5.5		2018-04-02 Final				

Labs, Investigations and Consults – no changes from current NL CPP and Other Sidebars.

Visits (Most Recent 10)

Updated to limit to most recent 10 visits and to include the **Concern** (if included on Appointment and/or Visit) and **Reason for Visit** (Diagnosis, if included Visit Note) and **Attending Provider**:



Immunizations – <u>no change</u> from current NL CPP and Other Sidebars.

Imported Hx

Will display any Task with a Category of Imported Hx and will display the scanned attachment icon, if Task includes a scanned attachment.

Imported Hx	≡
Complete Results Description	
Imported Hx, Test scanned CPP	<u>></u>